**cercon<sup>®</sup> TCT**

## **Cercon<sup>®</sup> TCT Stains**

Ergänzung zur Gebrauchsanweisung  
Directions for Use – Supplement  
Complément au mode d'emploi  
Integrazione delle istruzioni per l'uso  
Ampliación de las instrucciones de uso

**DENTSPLY**

# Cercon® TCT Stains

|                                         |    |
|-----------------------------------------|----|
| Ergänzung zur Gebrauchsanweisung        | 4  |
| Directions for Use – Supplement         | 7  |
| Complément au mode d’emploi             | 10 |
| Integrazione delle istruzioni per l’uso | 13 |
| Ampliación de las instrucciones de uso  | 16 |

# **Cercon® TCT Malfarben**

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## **Directions for Use – Supplement**

### **Indications:**

- Cercon TCT Stains and Glaze are indicated for staining and glazing Cercon ht restorations. They can also be used for staining and glazing with most lithium silicate based restorations, zirconia based restorations and most of the DENTSPLY Porcelains.

### **Contraindications:**

Only the indications listed above are suitable.

### **Warnings/Adverse Reactions:**

Consult Cercon TCT Stains and Glaze Safety Data Sheets (SDS) for safe handling and usage. Use only according to the indications listed above.

### **Precautions:**

Refer to the special notes throughout the instructions below for best practices for using Cercon TCT Stains and Glaze.

- This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.

### **Storage conditions:**

Store materials in packages with the lid tightly sealed. Do not leave open bottles where dust and debris may contaminate the porcelain. Replace cap when not in use. Avoid vibration and heat. Keep dry. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence. Do not use after expiration date.

# Firing chart

Please observe the instructions of our zirconium oxide discs Cercon® ht. There you will find information on the intended use, contraindications, warnings, cautions and technical data.

Make sure to apply the stains in very thin layers. For darker shades, an extra firing step may be added if necessary in order to obtain even better shade fidelity.

| °C                               | Dry   | Preheat | Low Temp | Heat Rate | High Temp | Vac Start Temp | Vac Stop Temp | Vakuum Zeit | Hold Time Air | Cool Time |
|----------------------------------|-------|---------|----------|-----------|-----------|----------------|---------------|-------------|---------------|-----------|
|                                  | (min) | (min)   | °C       | °C        | °C        | (min)          | (min)         | (min)       | (min)         | (min)     |
| Stains/Glaze → PFZ               | 3     | 3       | 450      | 55        | 820       | 0              | 0             | 0           | 1.5           | 6         |
| Stains/Glaze → Cercon ceram Kiss | 3     | 3       | 450      | 55        | 800       | 0              | 0             | 0           | 1.5           | 6         |
| Stains/Glaze → Cercon ht         | 3     | 3       | 450      | 55        | 820       | 0              | 0             | 0           | 1.5           | 6         |

| °F                               | Dry   | Preheat | Low Temp | Heat Rate | High Temp | Vac Start Temp | Vac Stop Temp | Hold Time | Hold Time Air | Cool Time |
|----------------------------------|-------|---------|----------|-----------|-----------|----------------|---------------|-----------|---------------|-----------|
|                                  | (min) | (min)   | °F       | °F        | °F        | (min)          | (min)         | (min)     | (min)         | (min)     |
| Stains/Glaze → PFZ               | 3     | 3       | 842      | 99        | 1,508     | 0              | 0             | 0         | 1.5           | 6         |
| Stains/Glaze → Cercon ceram Kiss | 3     | 3       | 842      | 99        | 1,472     | 0              | 0             | 0         | 1.5           | 6         |
| Stains/Glaze → Cercon ht         | 3     | 3       | 842      | 99        | 1,508     | 0              | 0             | 0         | 1.5           | 6         |

| Cercon TCT Stain & Glaze                            |       |  | REF     |                           |    |  | REF     |
|-----------------------------------------------------|-------|--|---------|---------------------------|----|--|---------|
| Cercon TCT Stain Kit                                |       |  |         | Cercon TCT Stain 4        | 5g |  | D651504 |
| 1x Universal Stain & Glaze Liquid                   |       |  | D651590 | Cercon TCT Stain i1       | 5g |  | D651511 |
| 1x Cercon TCT Glaze                                 |       |  |         | Cercon TCT Stain i2       | 5g |  | D651512 |
| 15 x Cercon TCT Stains                              |       |  |         | Cercon TCT Stain white    | 5g |  | D651520 |
|                                                     |       |  |         | Cercon TCT Stain creme    | 5g |  | D651521 |
| DENTSPLY Prosthetics Universal Stain & Glaze Liquid | 15 ml |  | D601315 | Cercon TCT Stain sunset   | 5g |  | D651522 |
| DENTSPLY Prosthetics Universal Stain & Glaze Liquid | 50 ml |  | D601350 | Cercon TCT Stain copper   | 5g |  | D651523 |
| Cercon TCT Glaze                                    | 5g    |  | D651322 | Cercon TCT Stain khaki    | 5g |  | D651524 |
| Cercon TCT Stain 0                                  | 5g    |  | D651500 | Cercon TCT Stain olive    | 5g |  | D651525 |
| Cercon TCT Stain 1                                  | 5g    |  | D651501 | Cercon TCT Stain mahogany | 5g |  | D651526 |
| Cercon TCT Stain 2                                  | 5g    |  | D651502 | Cercon TCT Stain purple   | 5g |  | D651505 |
| Cercon TCT Stain 3                                  | 5g    |  | D651503 |                           |    |  |         |

# The sleek color scheme

\* Cercon TCT Stains are necessary for individualisation not for coloring.

| V-SHADE | DISK-SHADE | 1ST STAIN FIRING                  | 2ND STAIN FIRING     |                       |
|---------|------------|-----------------------------------|----------------------|-----------------------|
|         |            | Stain/glaze                       | Stain/glaze          | Incisal               |
| A1      | A1*        | Stain 3, some olive               |                      | Stains i1 + i2        |
| A2      | A1         | Stains 2 + 3, some olive          |                      | Stains i1 + i2, white |
| A3      | A3*        | Stains 1 + 3                      |                      | Stains i1 + i2        |
| A3,5    | A3         | Olive                             | Stain 3              | Stains i1 + i2        |
| A4      | A3         | Stains 1 + 3                      | Stain 3              | Stains i1 + i2        |
| B1      | B1*        | Stain 3                           |                      | Stains i1 + i2        |
| B2      | B1         | Stains 3 +4                       |                      | Stain i1 + white      |
| B3      | B3         | Stains 2 + 3                      |                      | Stains i1 + i2        |
| B4      | B3*        | Stains 3 + 4, some olive + sunset |                      | Stains i1 + i2        |
| C1      | C1*        | Stain 3                           |                      | Stains i1 + i2        |
| C2      | C1         | Stain 3                           |                      | Stains i1 + i2        |
| C3      | C1         | Stain 3                           |                      | Stains i1 + i2        |
| C4      | C1         | Stain 3, some olive               |                      | Stains i1 + i2        |
| D2      | D2*        | Stain 4                           |                      | Stains i1 + i2        |
| D3      | D2         | Stain 4                           |                      | Stains i1 + i2        |
| D4      | D2         | Stain 3, some olive               | Stain 3, some sunset | Stains i1 + i2        |

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DeguDent GmbH
DENTSPLY International Inc.